

# $\delta^{18}\text{O}$ isotopic analysis of inorganic and organic solids by Thermal Conversion

Isotope Science Laboratory  
Department of Physics and Astronomy  
University of Calgary

## Overview:

The  $^{18}\text{O}/^{16}\text{O}$  ratio of pure inorganic and organic materials (e.g.  $\text{BaSO}_4$  or plant cellulose) is determined using a high temperature reactor coupled to a Continuous Flow - Isotope Ratio Mass Spectrometer (CF-IRMS). The system is comprised of a HEKA High Temperature Oxygen analyzer (HTO)<sup>®</sup> interfaced to a Thermo Delta+XL<sup>®</sup> mass spectrometer via a Conflow-III<sup>®</sup> open split/interface. Oxygen-containing compounds are quantitatively converted to CO at temperatures between 1350°C and 1450°C. After GC separation, the CO gas is swept by a carrier gas into the mass spectrometer and raw  $\delta^{18}\text{O}$  values are determined. Results are expressed in the usual per mil notation relative to the international V-SMOW standard. The instruments are fully automated and computer controlled using ISODAT 2.0 software. Prior to sample analysis all instrument conditions (reactor temperatures, GC oven temperature, carrier gas flow, ion beam background(s), ion source stability and signal linearity) are checked and recorded ("*Daily check*" routine).

Stable isotope ratios are expressed as delta ( $\delta$ ) and are measures of a 'per mill' (‰), or parts per thousand difference between the isotope ratio of a sample and that of a known (International) standard material.

## Online (CF-IRMS) technique:

1. Approximately 300µg of pure, homogenized sample are weighed into high purity silver cups (Elemental Microanalysis p/n: D2000)
2. The sample-containing cups are dropped into the high temperature (HTO) pyrolysis reactor using a Costech Zero Blank<sup>®</sup> autosampler
3. The HTO is maintained at temperatures between 1350°C and 1450°C for organic (e.g. cellulose) and inorganic (e.g.  $\text{BaSO}_4$ ) materials respectively
4. Carrier flow is ~90 ml/min of UHP 5.0 helium
5. The carrier gas sweeps the gaseous pyrolysis products through a GC column (Elemental MicroAnalysis p/n: E3042) in which CO is separated from  $\text{N}_2$
6. CO is introduced via an open split/interface directly into the ion source of the mass-spectrometer

7. Ion currents of masses 28 and 30 are measured simultaneously and the 30/28 ratio of the sample gas is compared to that of a standard CO reference (GTS-Praxair p/n: CO 4.7Z-AS)
8. Internal lab standards are used at the beginning, between (~ every 5th) and the end of each sequence to correct for instrument drift and to normalize the data to the international VSMOW-SLAP scale
9. Six (6) replicates of different weights are placed near the beginning of each sequence to allow for element wt.% determination and 'non-linearity' correction
10. Internal lab standards have been characterized against the International Standards listed below and are re-checked periodically
11. USGS LIMS is used for drift correction, normalization and data management (<http://water.usgs.gov/software/LIMS/>)

International Reference Materials for  $^{18}\text{O}$  solids:

<u>Material</u>	<u><math>\delta^{18}\text{O}_{\text{VSMOW}}</math></u>
IAEA NO <sub>3</sub>	+25.32 ±0.29‰
USGS34	-27.78 ±0.37‰
USGS35	+56.81 ±0.31‰
NBS 127	+8.59 ±0.26‰
IAEA SO5	+12.13 ±0.33‰
IAEA SO6	-11.35 ±0.31‰
IAEA 601	+23.14 ±0.19‰
IAEA 602	+71.28 ±0.36‰
VSMOW2	0‰
SLAP	-55.5‰

Accuracy and Precision:

Accuracy and precision of  $\delta^{18}\text{O}$  determinations on cellulose, AgNO<sub>3</sub> and BaSO<sub>4</sub> are generally better than ±0.3‰ (one standard deviation based on n=50 lab standards).

## References:

1. R. A. Werner, B. E. Kornexl, A. Roßmann, H. L. Schmidt., On-line determination of delta O-18 values of organic substances. *Anal. Chim. Acta* 1996, 319,159
2. G. D. Farquhar, K. B. Henry, J. M. Styles. A rapid on-line technique for determination of oxygen isotope composition of nitrogen-containing organic matter and water. *Rapid Commun. Mass Spectrom.* 1997, 11, 155
3. J. Koziat. Isotope ratio mass spectrometric method for the on-line determination of oxygen-18 in organic matter. *J. Mass Spectrom.* 1997, 32, 103
4. B.E. Kornexl, M. Gehre, R. Hofling, R. A. Werner. On-line  $\delta^{18}\text{O}$  measurement of organic and inorganic substances. *Rapid Communications in Mass Spectrometry* 1999, 13, 1685
5. M. Gehre, G. Strauch. High-temperature elemental analysis and pyrolysis techniques for stable isotope analysis. *Rapid Commun. Mass Spectrom.* 2003, 17, 1497
6. R.A. Werner and W.A. Brand, Referencing strategies and techniques in stable isotope ratio analysis, *Rapid Communications in Mass Spectrometry*, 2001; 15: 501-519
7. Böhlke J.K., Mroczkowski S.J. and Coplen T.B., Oxygen isotopes in nitrate: New reference materials for  $^{18}\text{O}:^{17}\text{O}:^{16}\text{O}$  measurements and observations on nitrate - water equilibration. *Rapid Communications in Mass Spectrometry*, 2003; 17: 1835-1846
8. W.E. Brand et al., Comprehensive inter-laboratory calibration of reference materials for  $\delta^{18}\text{O}$  versus VSMOW using various on-line high-temperature conversion techniques. *Rapid Communications in Mass Spectrometry*, 2009; 23: 999-1019
9. F. Accoe, M. Berglund, B. Geypens, P. Taylor. Methods to reduce interference effects in thermal conversion elemental analyzer/continuous flow isotope ratio mass spectrometry  $\delta^{18}\text{O}$  measurements for nitrogen-containing compounds. *Rapid Commun. Mass Spectrom.* 2008, 22, 2280
10. H. P. Sieper, H. J. Kupka, L. Lange, A. Roßmann, N.Tanz, H. L. Schmidt. Essential methodological improvements in the oxygen isotope ratio analysis of N-containing organic compounds. *Rapid Commun. Mass Spectrom.* 2010, 24, 2849
11. H. P. Qi, M. Gröning, T. B. Coplen, B. Buck, S. J. Mroczkowski, W. A. Brand, H. Geilmann , M. Gehre. Novel silver tubing method for quantitative introduction of water into high temperature conversion systems for stable hydrogen and oxygen isotopic measurements. *Rapid Commun. Mass Spectrom.* 2010, 24, 1821
12. Qi, H., Coplen, T.B. & Wassenaar, L.I., Improved online  $\delta^{18}\text{O}$  measurements of nitrogen- and sulfur-bearing organic materials and a proposed analytical protocol. *Rapid Communications in Mass Spectrometry*, 2011; 25: 2049– 2058
13. G.B. Hunsinger and L.A. Stern, Improved accuracy in high-temperature conversion elemental analyzer  $\delta^{18}\text{O}$  measurements of nitrogen-rich organics. *Rapid Commun. Mass Spectrom.* 2012, 26, 554–562
14. Roland A. Werner, The online  $^{18}\text{O}/^{16}\text{O}$  analysis: Development and application, *Isotopes Environ. Health Stud.*, Vol. 39, No. 2, June 2003, pp. 85–104