

Rotary Evaporator

Apparatus: IKA RV 10 Rotary Evaporator

Description: The rotovap can be used to evaporate solvent from a solid or liquid solute at a decreased temperature by heating with a universal water heating bath and applying a vacuum.

Specifications:

- Maximum temperature of 180 °C
- Maximum volume of 1L per sample

Rates (per sample):

\$10 (internal)

\$20 (external, academic)

\$30 (Industrial)

Contact information:

Hao Xu

Email: hao.xu1@ucalgary.ca

Lab: ENG 403